

Lecture #27: Quantum Dots (DQ) Imaging for Thermal flow studies

Dr. Hui HU

Martin C. Jischke Professor and Director

Advanced Flow Diagnostics and Experimental Aerodynamics Laboratory

Department of Aerospace Engineering, Iowa State University

2251 Howe Hall, Ames, IA 50011-2271

Email: huhui@iastate.edu

❑ Introduction about Quantum Dots


Quantum dots are *chemically synthesized semiconductor nanoparticles* that can exhibit *unique quantum confinement behavior*.

- *Size dependent emission:*
 - Smaller particles emit high-energy (shorter wavelength) radiation
 - Larger particles emit low-energy (longer wavelength) radiation.
- *Adjustable hydrophilic or hydrophobic properties:*
 - The surface layer can be designed to be either hydrophilic or hydrophobic allowing the nanoparticles to be solubilized in a variety of liquids (i.e. polar and nonpolar).



*The structure of a
(CdSe)ZnS quantum dot*

What's Quantum Dots




2.3nm → 5.5nm

□ ABSORPTION AND EMISSION SPECTRA OF (CdSe)ZnS QUANTUM DOTS



- **Absorption spectrum:**
 - *wide range of absorption*
 - *various laser can be used as excitation source*
- **Emission spectrum:**
 - *emission peak wavelength is adjustable by controlling the size of QDs*
 - *emission spectrum is symmetrical to the emission peak*
 - *emission concentrates near emission peak with FWHM only 30~40 nm*

Quantum Dots	Emission	
	Emission Peak	Full-Width Half-Maximum
QD-1	507 nm	33 nm
QD-2	593 nm	34 nm
QD-3	630 nm	36 nm
QD-4	650 nm	40 nm






❏ QUANTUM DOTS (QDs) FOR THERMO-FLUIDS MEASUREMENTS

- The propose of the present study is **to explore the possibility** of Quantum dots (QDs) for thermo-fluids diagnostics.
- Compared to commonly-used fluorescent dyes such as Fluorescein and Rhodamine B, **the photophysical properties** of (CdSe)ZnS QDs will be investigated:
 - Absorption and emission spectra.
 - Stability against photobleaching under different laser excitations.
 - Temperature sensitivity under different laser excitations.
- The **applications** of QDs for thermal-fluids diagnostics
 - **Individual QDs as particle tracers for velocimetry in micro flows:**
 - **Single Quantum Dot (QD) Velocimetry by Pouya et al. Experiments in Fluids, 2006.**
 - QDs as **fluorescent tracers:**
 - **Flow visualization and concentration measurements in a pulsed jet flow.**
 - **Temperature distribution mapping in a stratified flow.**



□ QD imaging for flow velocity measurements



- QD is much brighter than fluorescent nanoparticles for Micro-PIV measurements

- Image of 6 nm (CdSe)ZnS quantum dots in aqueous solution within 100 nm of the surface. Field of view is 153 μm·93 μm







Fig. 4 Histogram of measured streamwise (U -component) and cross-stream (V -component) velocities in the illuminated plane

- Shahram Pouya, Manoochehr Koochesfahani Preston Snee, Mounji Bawendi, Daniel Nocera, "Single quantum dot (QD) imaging of fluid flow near surfaces", *Experiments in Fluids* (2005) 39: 784–786

☐ PHOTOPHYSICAL STABILITY AGAINST PHOTBLEACHING



- *Same volume (200ml) of QD, Fluorescein and Rhodamine B solutions were prepared.*
- *Concentration of the Fluorescein and Rhodamine B aqueous solution is 1.0×10^{-6} M.*
- *The solvent of QD is Hexane.*
- *The absolute molar concentration of QD is unknown, but it has the same initial fluorescence emission intensity as the Rhodamine B - 1.0×10^{-6} M solution.*



Experimental setup for the photoluminescence stability test



PHOTOLUMINESCENCE STABILITY AGAINST PHOTOBLEACHING (COMPARISON OF QD WITH FLUORESCEIN AND RHODAMINE B)

excitation laser source: pulsed excimer UV laser (308nm)




PHOTOLUMINESCENCE STABILITY AGAINST PHOTOBLEACHING (COMPARISON OF QD WITH FLUORESCEIN AND RHODAMINE B)

excitation laser source: Argon-ion laser (514.5 nm)




Temperature Sensitivity Test




Experimental setup for temperature sensitivity test

TEMPERATURE SENSITIVE OF QD, FLUORESC EIN AND RHODAMINE B


QD



Rhodamine B




Fluorescein



	308nm Excimer laser	488nm Argon-ion laser	514nm Argon-ion laser
QD-3	-0.94% per °C (5°C ~ 65°C)	-1.14% per °C (5°C ~ 65°C)	-1.05% per °C (5°C ~ 65°C)
Rhodamine B	-1.57% per °C (5°C ~ 65°C)	-1.20% per °C (5°C ~ 65°C) Kim & Khim (2001) -1.35 % per °C (15°C ~ 40°C)	-1.52% per °C (5°C ~ 65°C) Coppeta & Roger(1998) -1.54 % per °C (20°C ~ 60°C)
Fluorescein	-0.45% per °C (5°C ~ 65°C)	-0.25% per °C (5°C ~ 65°C) Coppeta & Roger(1998) -0.16 % per °C (20°C ~ 60°C)	+2.25% per °C (5°C ~ 65°C) Coppeta & Roger(1998) +2.43 % per °C (20°C ~ 60°C)

APPLICATION QUANTUM DOTS FOR THERMAL-FLUIDS DIAGNOSTICS (CONCENTRATION MEASUREMENTS IN A PULSED JET)



APPLICATION QUANTUM DOTS FOR THERMAL-FLUIDS DIAGNOSTICS (CONCENTRATION MEASUREMENTS IN A PULSED JET)


concentration

0.0


1.0



a. $t = t_0$




b. $t = t_0 + 0.2 \text{ s}$



c. $t = t_0 + 0.4 \text{ s}$







APPLICATION QUANTUM DOTS FOR THERMAL-FLUIDS DIAGNOSTICS

(TEMPERATURE MEASUREMENTS IN A STRATIFIED FLOW)





APPLICATION QUANTUM DOTS FOR THERMAL-FLUIDS DIAGNOSTICS

(TEMPERATURE MEASUREMENTS IN A STRATIFIED FLOW)



APPLICATION QUANTUM DOTS FOR THERMAL-FLUIDS DIAGNOSTICS

(TEMPERATURE MEASUREMENTS IN A STRATIFIED FLOW)



Temperature profiles along the vertical line in the center of the test cell.